

CERTIFICATE OF ANALYSIS

Anaplasma phagocytophilum

FA Positive Control

Catalog No.:	PC-IFA-AP
Lot:	P170517-001
Expiration:	26 May 2021

Quality Control Method: Indirect FA using Anaplasma phagocytophilum 12-well slide (catalog no.

SLD-IFA-AP), negative control (catalog no. NC-IFA-AP), anti-equine IgG

FITC conjugate (catalog no. CJ-F-EQUG-1ML or 10ML).

Specific Reaction: 2-4+ signal with the positive control at neat, no background. Trace-1+

endpoint reached at a 1/8 dilution.

Other Comments: NA

Will Thurman QC Manager:

QC Release Date: 26 May 2017

Pattern of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered

elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: Version #1