

CERTIFICATE OF ANALYSIS

Anaplasma phagocytophilum

FA Positive Control

Catalog No.:	PC-IFA-AP
Lot:	P191230-001
Expiration:	04 June 2023

Quality Control Method: Indirect FA using *Anaplasma phagocytophilum* 12-well slide (catalog no. SLD-IFA-AP), negative control (catalog no. NC-IFA-AP), anti-equine IgG FITC conjugate (catalog no. CJ-F-EQUG-1ML or 10ML).

Specific Reaction: 3-4+ signal with the control at neat, no background. Endpoint is 1-2+ at 1/8 dilution.

Other Comments: NA

QC Manager: David Hyatt



QC Release Date: 04 June 2019

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

27 June 2019