

CERTIFICATE OF ANALYSIS

Anaplasma phagocytophilum

FA Positive Control

Catalog No.:	PC-IFA-AP
Lot:	P220728-001
Expiration:	11 August 2030

Quality Control Method: Indirect FA using *Anaplasma phagocytophilum* 12-well slide (catalog no. SLD-IFA-AP), negative control (catalog no. NC-IFA-AP), and anti-equine IgG FITC conjugate (catalog no. CJ-F-EQUG-10ML).

Specific Reaction: 2–4+ signal with the control at neat. No background. Endpoint signal was reached at 1/8 dilution.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

19 August 2022