

**CERTIFICATE OF ANALYSIS**

*Brucella canis*

FA Positive Control

Catalog No.:	PC-IFA-CB
Lot:	P150813-001
Expiration:	01 September 2019

Quality Control Method: IFA using *Brucella canis* 12-well slide (catalog no. SLD-IFA-CB), negative control (catalog no. NC-IFA-CB), and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 3+ signal with the positive control, negative with the negative control. 1+ endpoint reached with 3 two-fold dilutions.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 01 September 2015

Pattern Of Fluorescence: Brightly staining short bacilli against uniform background of ovine red blood cells.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COAs and TDS\COA (VMRD)\Positive Controls\Canine Brucellosis\PC-IFA-CB\_P150813001\_190901.docx

02 September 2015