

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV) IgG

FA Positive Control

Catalog No.:	PC-IFA-CDV-G
Lot:	P230127-001
Expiration:	02 February 2031

Quality Control Method: Indirect FA on CDV 12-well slide (catalog no. SLD-IFA-CDV), and anti-canine IgG affinity purified FITC conjugates (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 3-4+ signal with the control at neat, no background. Endpoint signal was greater than 1/16 dilution.

Other Comments: NA

QC Manager: 

Jennifer Boren

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane fluorescence. High titer test sera will fluoresce brighter than this positive control serum.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

08 February 2023