

CERTIFICATE OF ANALYSIS

Feline Calicivirus (FCV)


FA Positive Control

Catalog No.:	PC-IFA-FCV
Lot:	P191004-003
Expiration:	04 April 2031

Quality Control Method: Indirect FA using FCV 12-well slide (catalog no. SLD-IFA-FCV), negative control (catalog no. NC-IFA-FCV), and anti-feline IgG affinity purified conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: 3–4+ signal with the positive control at neat. No background. Endpoint is at a dilution greater than 1/16.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Diffuse and granular cytoplasmic fluorescence with brightly fluorescing nuclear membranes, also brightly fluorescing individual nuclei.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

6 April 2023