

CERTIFICATE OF ANALYSIS

Feline Calicivirus (FCV)

FA Positive Control

| | |
|--------------|--------------|
| Catalog No.: | PC-IFA-FCV |
| Lot: | P230606-001 |
| Expiration: | 08 June 2031 |

Quality Control Method: Indirect FA using FCV 12-well slide (catalog no. SLD-IFA-FCV), negative control (catalog no. NC-IFA-FCV), and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: 2–4+ signal with the control at neat. No background. Endpoint is at a dilution greater than 1/16.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Diffuse and granular cytoplasmic fluorescence with brightly fluorescing nuclear membranes, also brightly fluorescing individual nuclei.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

12 June 2023