

CERTIFICATE OF ANALYSIS

Feline Infectious Peritonitis Virus Type 1 (FIP-1)


FA Positive Control

Catalog No.:	PC-IFA-FIP1
Lot:	P160606-002
Expiration:	16 April 2024

Quality Control Method: Indirect FA using FIP-1 12-well slide (catalog no. SLD-IFA-FIP1), negative control (catalog no. NC-IFA-FIP1), and anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: Trace-2+ signal with the positive control at neat with no background. 1+ endpoint is reached after 3 twofold dilutions.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern of Fluorescence: Multiple syncytia, plaques and some individual cells with granular cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Positive Controls\FIP-1\PC-IFA-FIP1_P160606-002_240416.docx

21 April 2020