

## **CERTIFICATE OF ANALYSIS**

## Feline Infectious Peritonitis Virus Type 1 (FIP-1)

FA Positive Control

Catalog No.:	PC-IFA-FIP1
Lot:	P210810-001
Expiration:	20 May 2025

Quality Control Method: Indirect FA using FIP-1 12-well slide (catalog no. SLD-IFA-FIP1), negative

control (catalog no. NC-IFA-FIP1), and anti-feline IgG FITC affinity purified

conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: 3-4+ signal with the control at neat, no background. Endpoint is greater

than four 2-fold dilutions.

Other Comments: NA

QC Manager:

David Hyatt

Havel Wyath

Pattern Of Fluorescence: Multiple syncytia, plaques and some individual cells with granular

cytoplasmic fluorescence.

Technical Data Sheet Version: 2