

CERTIFICATE OF ANALYSIS

Leishmania infantum

FA Positive Control

Catalog No.:	PC-IFA-LSH
Lot:	P160302-002
Expiration:	07 April 2024

Quality Control Method: Indirect FA using *Leishmania infantum* 12-well slide (catalog no. SLD-IFA-LSH), negative control (catalog no. NC-IFA-LSH), and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 1-3+ signal with the positive control at neat, trace background. Trace-2+ endpoint reached in 4 twofold dilutions.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Diffuse membranous fluorescence of organisms varying in shape from globus to trypanosome-like, some with clearly defined flagella. Single organisms and large clumps. High titer test sera will fluoresce brighter than this positive control serum.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\Positive Controls\Leishmania\PC-IFA-LSH_P160302-002_240407.docx

15 April 2020