

CERTIFICATE OF ANALYSIS

Leishmania infantum

FA Positive Control

Catalog No.:	PC-IFA-LSH
Lot:	P200319-001
Expiration:	06 March 2032

Quality Control Method: Indirect FA using Leishmania infantum 12-well slide (catalog no. SLD-IFA-

LSH), negative control (catalog no. NC-IFA-LSH), and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 2-4+ signal with the control at neat. No background. Endpoint dilution is

greater than 1/16.

Other Comments: NA

QC Manager:

David Hyatt

Waved Wyath

Pattern Of Fluorescence: Diffuse membranous fluorescence of organisms varying in shape from

globus to trypanosome-like, some with clearly defined flagella. Single organisms and large clumps. High titer test sera will fluoresce brighter than

this positive control serum.

Technical Data Sheet Version: 1

FOR IN VITRO LABORATORY USE ONLY.