

CERTIFICATE OF ANALYSIS

*Neospora caninum*


FA Positive Control

Catalog No.:	PC-IFA-NC-BOV
Lot:	P230228-001
Expiration:	09 March 2031

**Quality Control Method:** Indirect FA using *Neospora caninum* 12-well slide (catalog no. SLD-IFA-NC), bovine negative control (catalog no. NC-IFA-NC-BOV), and anti-bovine IgG1, 2 FITC affinity purified conjugate (catalog no. CJ-F-BOVG-AP-1ML or 10ML).

**Specific Reaction:** 3-4+ signal with the control at neat, no background. Endpoint signal was reached at 1/16 dilution.

**Other Comments:** When using *Neospora caninum* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

**QC Manager:**   
Jennifer Boren

**Pattern Of Fluorescence:** Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

16 March 2023