

CERTIFICATE OF ANALYSIS

Neospora caninum


FA Positive Control

Catalog No.:	PC-IFA-NC-CAN
Lot:	P210721-001
Expiration:	09 September 2025

Quality Control Method: Indirect FA using *Neospora caninum* 12-well slide (catalog no. SLD-IFA-NC), canine negative control (catalog no. NC-IFA-NC-CAN), and anti-canine FITC AP conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 2–4+ signal with the control at neat, no background. Endpoint signal was reached at 1/16 dilution.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: 1