

CERTIFICATE OF ANALYSIS

Neospora caninum


FA Positive Control

Catalog No.:	PC-IFA-NC-CAN
Lot:	P220601-001
Expiration:	14 June 2031

Quality Control Method: Indirect FA using *Neospora caninum* 12-well slide (catalog no. SLD-IFA-NC), canine negative control (catalog no. NC-IFA-NC-CAN), and anti-canine affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 4+ signal with the control at neat. No background. Endpoint is at a dilution greater than 1/16.

Other Comments: When using *Neospora caninum* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: 1