

CERTIFICATE OF ANALYSIS

Porcine Circovirus Type 2 (PCV-2)

FA Positive Control

Catalog No.:	PC-IFA-PCV2
Lot:	P150506-001
Expiration:	16 March 2031

Quality Control Method: Indirect FA using PCV-2 12-well slide (catalog no. SLD-IFA-PCV2), negative control (catalog no. NC-IFA-PCV2), and anti-porcine IgG affinity purified conjugate (catalog no. CJ-F-PORG-AP-10ML).

Specific Reaction: 3-4+ signal with the control at neat, no background. Endpoint signal was reached at 1/16 dilution.

Other Comments: NA

QC Manager: 
Jennifer Boren

Pattern Of Fluorescence: Intense whole nuclear and variable sized but discrete cytoplasmic inclusions.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

16 March 2023