

CERTIFICATE OF ANALYSIS

Porcine Epidemic Diarrhea Virus (PEDV)


FA Positive Control

Catalog No.:	PC-IFA-PEDV
Lot:	P141119-001
Expiration:	31 August 2030

Quality Control Method: Indirect FA using PEDV substrate slide (catalog no. SLD-IFA-PEDV), negative control (catalog no. NC-IFA-PEDV), and anti-porcine IgG affinity purified conjugate (catalog no. CJ-F-PORG-AP-1ML or 10ML).

Specific Reaction: 2–3+ signal at neat, trace background. An endpoint titer of trace to 1+ was reached at 1/8 dilution.

Other Comments: N/A

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Brightly fluorescing irregular cytoplasmic inclusions, mainly in single cells.

Technical Data Sheet Version: 1