

CERTIFICATE OF ANALYSIS

Toxoplasma gondii IgG

FA Positive Control

Catalog No.:	PC-IFA-TOXO-FEL-G
Lot:	P180724-001
Expiration:	14 August 2022

Quality Control Method: Indirect FA using *Toxoplasma gondii* 12-well slide (catalog no. SLD-IFA-TOXO), negative control (catalog no. NC-IFA-TOXO-FEL), and anti-feline IgG FITC conjugate (catalog no. CJ-F-FELG-1ML or 10ML).

Specific Reaction: 2 – 4+ signal with the control at neat, no background. Endpoint is trace – 1+ at 1/16 dilution.

Other Comments: When using *Toxoplasma gondii* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

QC Manager: Chris Cay



QC Release Date: 17 August 2018

Pattern Of Fluorescence: Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

24 August 2018