

CERTIFICATE OF ANALYSIS

Toxoplasma gondii


FA Positive Control

Catalog No.:	PC-IFA-TOXO-FEL-G
Lot:	P220606-001
Expiration:	15 June 2026

Quality Control Method: Indirect FA using *Toxoplasma gondii* 12-well slide (catalog no. SLD-IFA-TOXO), negative control (catalog no. NC-IFA-TOXO-FEL), and anti-feline IgG FITC conjugate (catalog no. CJ-F-FELG-10ML).

Specific Reaction: 3–4+ signal with the control at neat. No background. Endpoint signal of 1–2+ was reached at 1/16 dilution.

Other Comments: When using *Toxoplasma gondii* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

29 June 2022