



**VMRD, Inc.**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: [vtech@vmrd.com](mailto:vtech@vmrd.com)

Web site: <http://www.vmrd.com>

---

## Certificate of Analysis

### **POWDERED FA RINSE BUFFER**

**CATALOG NO.:** 210-90-RB

**VOLUME:** 1 pouch

**LOT:** P031022-001

**EXPIRATION:** 7 November 2005

**DESCRIPTION:** FA Rinse Buffer Powder, pH 9.0-9.5

**INTENDED USE:** FA Rinse Buffer is used to rinse FA substrate slides.

**RECONSTITUTION AND STORAGE:** Dissolve in 1 liter of deionized/distilled water to make 4X concentrate. Dilute this concentrate 1:4 for final use. Powdered and reconstituted FA Rinse Buffer should be stored at room temperature.

FOR RESEARCH AND INVESTIGATIONAL USE ONLY.
--

**WARRANTY:** *VMRD, Inc.* warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. *VMRD, Inc.* is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.