

Feline Leukemia Virus (FeLV)

Negative Blood Smear Slide

Catalog No.:	SLD-BSN-FELV
Shelf Life:	4 years from date of qualification
Substrate:	Blood negative for Feline Leukemia Virus (FeLV)
Species of Origin:	Feline

Description:

Feline blood smeared on glass slide negative for FeLV.

Intended Use:

As a fluorescent antigen negative control slide.

Storage:

Store sealed in foil pouch at $<-10^{\circ}\text{C}$. Avoid self-defrosting freezers.

References:

NA

Technical Data Sheet Version:

Version 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\TDS\Blood Smears\SLD-BSN-FELV TDS V2 CPM Active MKT.DOCX

04 April 2017

