

CERTIFICATE OF ANALYSIS

Bovine Adenovirus Type 3 (BAV-3)

FA Control Slide

Catalog No.:	SLD-FAC-BAV3
Lot:	P161010-001
Expiration:	30 November 2024

Quality Control Method: Direct FA using BAV-3 FITC conjugate (catalog no. CJ-F-BAV3-1ML or 10ML).

Specific Reaction: 2-4+ signal with no background on the positive well. Negative on the negative well. There are 0-50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern of Fluorescence: Single cells with small granular nuclear inclusions.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\2-well Slides\BAV3\SLD-FAC-BAV3_P161010-001_241130.docx

04 December 2020