

CERTIFICATE OF ANALYSIS

Bovine Adenovirus Type 3 (BAV3)

FA Control Slide

Catalog No.:	SLD-FAC-BAV3
Lot:	P210320-001
Expiration:	26 April 2025

Quality Control Method: Direct FA using BAV-3 FITC conjugate (catalog no. CJ-F-BAV3-10ML).

Specific Reaction: 2–4+ signal on the positive well and negative on the negative well with trace–1+ background. There are 0–25 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Single cells with small granular nuclear inclusions.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

05 May 2021