

CERTIFICATE OF ANALYSIS

Bovine Adenovirus Type 5 (BAV5)

FA Control Slide

Catalog No.:	SLD-FAC-BAV5
Lot:	P210325-001
Expiration:	24 June 2025

Quality Control Method: Direct FA using BAV-5 FITC conjugate (catalog no. CJ-F-BAV5-10ML).

Specific Reaction: 2-4+ signal on the positive well and negative on the negative well with no background. There are 0-25 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Diffuse and granular cytoplasmic fluorescence. Brightly fluorescing nuclear inclusions and nuclei. Some infected cells are multinucleated.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

30 June 2021