

CERTIFICATE OF ANALYSIS

Bovine Coronavirus (BCV)

FA Control Slide

Catalog No.:	SLD-FAC-BCV
Lot:	P150109-003
Expiration:	02 February 2019

Quality Control Method: Direct FA using BCV FITC conjugate (catalog no. CJ-F-BCV-1ML or 10ML).

Specific Reaction: 2-4+ fluorescence on the positive well and negative on the negative well with no background. There are 0 to too numerous to count infected cells per high-power field.

Other Comments: NA

QC Manager: *Torrey Beeson* Torrey Beeson

QC Release Date: 03 February 2015

Pattern Of Fluorescence: Granular cytoplasmic and nuclear fluorescence mainly on single cells but also some small syncytia.

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\COA (VMRD)\2-well Slides\BCV\SLD-FAC-BCV_P150109-003_190202.docx

04 February 2015