

CERTIFICATE OF ANALYSIS

Bovine Coronavirus (BCV)

FA Control Slide

Catalog No.:	SLD-FAC-BCV
Lot:	P171020-010
Expiration:	21 May 2025

Quality Control Method: Direct FA using BCV FITC conjugate (catalog no. CJ-F-BCV-1ML or 10ML).

Specific Reaction: 2-3+ signal on the positive well and negative on the negative well with trace background. There are 0 to 50 infected cells per high-power field.

Other Comments: This lot of slides exhibits inconsistent distribution of infected cells.

QC Manager:



David Hyatt

Pattern Of Fluorescence: Granular cytoplasmic and nuclear fluorescence mainly on single cells but also some small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\2-well Slides\BCV\SLD-FAC-BCV_P171020-010_250521.docx

01 June 2021