

CERTIFICATE OF ANALYSIS

## Bovine Coronavirus (BCV)

FA Control Slide

Catalog No.:	SLD-FAC-BCV
Lot:	P220412-001
Expiration:	21 November 2030

Quality Control Method: Direct FA using BCV FITC conjugate (catalog no. CJ-F-BCV-10ML).

Specific Reaction: 1–3+ signal on the positive well and negative on the negative well. Trace background. There are 0–50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Granular cytoplasmic and nuclear fluorescence mainly on single cells but also some small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

30 November 2022