

**CERTIFICATE OF ANALYSIS**

## Bovine Parvovirus (BPV)

FA Control Slide

Catalog No.:	SLD-FAC-BPV
Lot:	P161013-001
Expiration:	10 November 2020

Quality Control Method: Direct FA using BPV FITC conjugate (catalog no. CJ-F-BPV-1ML or 10ML).

Specific Reaction: 3-4+ signal on the positive well, negative on the negative well. No background. There are 0-25 infected cells per high-power field.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 10 November 2016

Pattern of Fluorescence: Ground glass and/or granular cytoplasmic fluorescence with occasional bright and diffuse nuclear fluorescence.

Technical Data Sheet Version: Version # 2

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Quality\COA (VMRD)\2-well Slides\BPV\SLD-FAC-BPV\_P161010-001\_201110.docx

11 November 2016