

CERTIFICATE OF ANALYSIS

# Bovine Respiratory Syncytial Virus (BRSV)

FA Control Slide

Catalog No.:	SLD-FAC-BRSV
Lot:	P200810-001
Expiration:	18 August 2025

Quality Control Method: Direct FA using BRSV FITC conjugate (catalog no. CJ-F-BRSV-10ML).

Specific Reaction: 2–4+ signal on the positive well and negative on the negative well with no background. There are 0–50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Granular cytoplasmic fluorescence with some small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

02 September 2021