

CERTIFICATE OF ANALYSIS

Bovine Viral Diarrhea Virus (BVDV)

FA Control Slide

Catalog No.:	SLD-FAC-BVD
Lot:	P181201-001
Expiration:	22 January 2023

Quality Control Method: Direct FA using BVDV FITC conjugate (catalog no. CJ-F-BVD-1ML or 10ML).

Specific Reaction: 2 – 3+ signal on the positive well and negative on the negative well with trace background. There are 1 - 50 infected cells per high-power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 24 January 2019

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

28 January 2019