

CERTIFICATE OF ANALYSIS

Bovine Viral Diarrhea Virus (BVDV)


FA Control Slide

Catalog No.:	SLD-FAC-BVD
Lot:	P190830-006
Expiration:	05 December 2023

Quality Control Method: Direct FA using BVDV FITC conjugate (catalog no. CJ-F-BVD-10ML).

Specific Reaction: 2-4+ signal with trace background on the positive well and negative with trace-1+ background on the negative well. There are 0-50 infected cells per high-power field.

Other Comments: NA

QC Manager: 
David Hyatt

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

8 January 2020