

CERTIFICATE OF ANALYSIS

Bovine Viral Diarrhea Virus (BVDV)

FA Control Slide

Catalog No.:	SLD-FAC-BVD
Lot:	P211205-002
Expiration:	29 December 2025

Quality Control Method: Direct FA using BVDV FITC conjugate (catalog no. CJ-F-BVD-10ML).

Specific Reaction: 2–4+ signal with trace background on the positive well and negative with trace background on the negative well. There are 0–50 infected cells per high-power field.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Individual cells with smooth, undifferentiated and/or “ground glass” cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

10 January 2022