

CERTIFICATE OF ANALYSIS

# Canine Coronavirus (CCV)

FA Control Slide

Catalog No.:	SLD-FAC-CCV
Lot:	P220702-001
Expiration:	02 September 2030

Quality Control Method: Direct FA using CCV conjugate (catalog no CJ-F-CCV-10ML).

Specific Reaction: 2-4+ signal on the positive well and negative on the negative well. Trace background. There are 0-20 infected cells per high-power field.

Other Comments: NA

QC Manager:   
David Hyatt

Pattern Of Fluorescence: Large fluorescent syncytia with some single cells in between.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

07 September 2022