

## **CERTIFICATE OF ANALYSIS**

## Canine Distemper Virus (CDV)

FA Control Slide

Catalog No.:	SLD-FAC-CDV
Lot:	P220315-001
Expiration:	27 April 2026

Quality Control Method: Direct FA on CDV polyclonal conjugate (catalog no. CJ-F-CDV-10ML) and

CDV monoclonal conjugate (catalog no. CJ-F-CDV-MAB-10ML).

Specific Reaction: 2–4+ signal on the positive well and negative on the negative well. No

background. There are 0-50 infected cells per high-power field.

Other Comments: NA

QC Manager:

David Hyatt

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane

fluorescence.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.