

CERTIFICATE OF ANALYSIS

Canine Parainfluenza Virus Type 2 (CPI-2)

FA Control Slide

Catalog No.:	SLD-FAC-CPI
Lot:	P220401-001
Expiration:	19 May 2026

Quality Control Method: Direct FA using CPI-2 direct conjugate (catalog no. CJ-F-CPI-10ML).

Specific Reaction: 3–4+ signal on the positive well. Negative with trace background on the negative well. There are 0–20 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Rough, fibrous and speckled cytoplasmic and undifferentiated nuclear fluorescence in syncytia and single cells.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

26 May 2022