

CERTIFICATE OF ANALYSIS

## Feline Calicivirus (FCV)

FA Control Slide


Catalog No.:	SLD-FAC-FCV
Lot:	P190405-001
Expiration:	21 March 2031

**Quality Control Method:** Indirect FA using FCV 2-well slide (catalog no. SLD-FAC-FCV), positive control (catalog no. PC-IFA-FCV), and anti-feline IgG affinity purified conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** 1–3+ signal on the positive well and negative on the negative well. No background. There are 0–50 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:** David Hyatt



**Pattern Of Fluorescence:** Diffuse and granular cytoplasmic fluorescence with brightly fluorescing nuclear membranes, also brightly fluorescing individual nuclei.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

31 March 2023