

CERTIFICATE OF ANALYSIS

# Feline Leukemia Virus (FeLV)


FA Control Slide

Catalog No.:	SLD-FAC-FELV
Lot:	P191104-001
Expiration:	08 August 2031

**Quality Control Method:** Indirect FA using primary anti-FeLV antibody (catalog no. AB1-FELV) and secondary anti-caprine IgG FITC conjugate (catalog no. AB2-FELV).

**Specific Reaction:** 2–4+ signal, trace–1+ background on the positive well and negative, trace–1+ background on the negative well. There are 0–25 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:**   
David Hyatt

**Pattern Of Fluorescence:** Individual cell cytoplasmic fluorescence with occasional small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

30 August 2023