

CERTIFICATE OF ANALYSIS

Feline Leukemia Virus (FeLV)


FA Control Slide

Catalog No.:	SLD-FAC-FELV
Lot:	P230123-001
Expiration:	13 March 2031

Quality Control Method: Indirect FA using primary anti-FeLV antibody (catalog no. AB1-FELV) and secondary anti-caprine IgG FITC conjugate (catalog no. AB2-FELV).

Specific Reaction: 3-4+ signal on the positive well and negative on the negative well with trace background. There are 0-10 infected cells per high-power field.

Other Comments: NA

QC Manager: 
Jennifer Boren

Pattern Of Fluorescence: Individual cell cytoplasmic fluorescence with occasional small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

14 March 2023