

CERTIFICATE OF ANALYSIS

Feline Panleukopenia Virus (FPLV)


FA Control Slide

Catalog No.:	SLD-FAC-FPL
Lot:	P190226-003
Expiration:	21 February 2031

Quality Control Method: Direct FA using FPLV FITC conjugate (catalog no. CJ-F-FPL-MAB-10ML).

Specific Reaction: 3–4+ signal on the positive well and negative on the negative well. Trace background. There are 20–60 infected cells per high-power field.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Single cells with cytoplasmic and/or nuclear fluorescence

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

28 February 2023