

CERTIFICATE OF ANALYSIS

Porcine Circovirus Type 2 (PCV-2)

FA Control Slide

Catalog No.:	SLD-FAC-PCV2
Lot:	P230403-001
Expiration:	18 April 2031

Quality Control Method: Direct FA using PCV-2 conjugate (catalog no. CJ-F-PCV2-10ML).

Specific Reaction: 3–4+ signal on the positive well and negative on the negative well. No background. There are 0–15 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Intense whole nuclear and variable sized but discrete cytoplasmic inclusions.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

20 April 2023