

CERTIFICATE OF ANALYSIS

Porcine Epidemic Diarrhea Virus (PEDV)


FA Control Slide

Catalog No.:	SLD-FAC-PEDV
Lot:	P150505-001
Expiration:	04 March 2032

Quality Control Method: Indirect FA using PEDV positive control (catalog no. PC-IFA-PEDV), and anti-porcine IgG affinity purified conjugate (catalog no. CJ-F-PORG-AP-10ML).

Specific Reaction: 2–3+ fluorescence on the positive well and negative on the negative well. No background. There are 0–25 infected cells per high-power field.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Brightly fluorescing irregular cytoplasmic inclusions, mainly in single cells.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

07 March 2024