

CERTIFICATE OF ANALYSIS

Parainfluenza Virus Type 3 (PI-3)

FA Control Slide

Catalog No.:	SLD-FAC-PI3
Lot:	P190429-001
Expiration:	15 May 2023

Quality Control Method: Direct FA using PI-3 conjugate (catalog no. CJ-F-PI3-1ML or 10ML).

Specific Reaction: 2 – 4+ signal on the positive well and negative on the negative well with no background. There are 0 – 30 infected per high-power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 15 May 2019

Pattern Of Fluorescence: Small, bright spherical cytoplasmic inclusions in syncytia and single cells.

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

24 May 2019