

CERTIFICATE OF ANALYSIS

Parainfluenza Virus Type 3 (PI-3)

FA Control Slide

| | |
|--------------|---------------|
| Catalog No.: | SLD-FAC-PI3 |
| Lot: | P220215-001 |
| Expiration: | 30 March 2026 |

Quality Control Method: Direct FA using PI-3 conjugate (catalog no. CJ-F-PI3-10ML).

Specific Reaction: 3–4+ signal on the positive well and negative on the negative well with no background. There are 0–10 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Small, bright spherical cytoplasmic inclusions in syncytia and single cells.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

05 April 2022