

CERTIFICATE OF ANALYSIS

# Porcine Parvovirus (PPV)

FA Control Slide

Catalog No.:	SLD-FAC-PPV
Lot:	P180130-001
Expiration:	02 February 2022

Quality Control Method: Direct FA using PPV conjugate (catalog no. CJ-F-PPV-1ML or 10ML).

Specific Reaction: 2-4+ signal with trace background on the positive well and negative on the negative well with no background. There are 0 - 50 infected cells per high-power field.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 06 February 2018

Pattern Of Fluorescence: Fluorescence limited almost entirely to nucleus of cell. Some granular cytoplasmic staining.

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\2-well Slides\PPV\SLD-FAC-PPV\_P180130-001\_220202.docx

20 February 2018