

CERTIFICATE OF ANALYSIS

# Porcine Parvovirus (PPV)

FA Control Slide

Catalog No.:	SLD-FAC-PPV
Lot:	P220215-003
Expiration:	31 March 2026

Quality Control Method: Direct FA using PPV conjugate (catalog no. CJ-F-PPV-10ML).

Specific Reaction: 1–3+ signal on the positive well and negative on the negative well. No background. There are 0–20 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Fluorescence limited almost entirely to nucleus of cell. Some granular cytoplasmic staining.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

04 April 2022