

CERTIFICATE OF ANALYSIS

Rabies Recombinant Nucleoprotein (rNP)

FA Control Slide

Catalog No.:	SLD-FAC-RAB
Lot:	P190422-001
Expiration:	15 March 2031

Quality Control Method: Direct FA using purchased Rabies FITC conjugate.

Specific Reaction: 2–4+ signal on the positive well and negative on the negative well. No background. There are 0–20 infected cells per high-power field.

Other Comments: NA

QC Manager: 
Jennifer Boren

Pattern Of Fluorescence: Smooth and particulate cytoplasmic fluorescence of single and multinucleated cells.

Technical Data Sheet Version: 3

FOR *IN VITRO* LABORATORY USE ONLY.

16 March 2023