

CERTIFICATE OF ANALYSIS

Rabies Recombinant Nucleoprotein (rNP)

FA Control Slide

Catalog No.:	SLD-FAC-RAB
Lot:	P190422-002
Expiration:	09 April 2023

Quality Control Method: Direct FA using Rabies FITC in-house conjugate.

Specific Reaction: 2-4+ signal on the positive well, negative on the negative well with no background. There are 0-20 infected cells per high powered field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 09 April 2019

Pattern of Fluorescence: Smooth and particulate cytoplasmic fluorescence of single and multinucleated cells.

Technical Data Sheet Version: Version #3

FOR *IN VITRO* LABORATORY USE ONLY.

16 April 2019