

CERTIFICATE OF ANALYSIS

Reovirus (REO)

FA Control Slide

Catalog No.:	SLD-FAC-REO
Lot:	P151013-001
Expiration:	22 May 2023

Quality Control Method: Direct FA using REO conjugate (catalog no. CJ-F-REO-1ML or 10ML).

Specific Reaction: 3 - 4+ signal on positive well, negative on negative well with no background noted on either well. 0-50 infected cells per high powered field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 22 May 2019

Pattern Of Fluorescence: Individual cells or plaques with diffuse cytoplasmic fluorescence and characteristic bright cytoplasmic inclusions ranging in size from medium to large.

Technical Data Sheet Version: Version #1

FOR *IN VITRO* LABORATORY USE ONLY.

24 May 2019