

CERTIFICATE OF ANALYSIS

# Transmissible Gastroenteritis Virus (TGEV)

FA Control Slide

Catalog No.:	SLD-FAC-TGE
Lot:	P220514-004
Expiration:	23 March 2026

Quality Control Method: Direct FA using TGEV conjugate (catalog no. CJ-F-TGE-10ML).

Specific Reaction: 3–4+ signal with no background on the positive well and negative with trace background on the negative well. There are 0–10 infected cells and/or syncytia per high-power field.

Other Comments: NA

QC Manager:   
David Hyatt

Pattern Of Fluorescence: Moderate to large syncytia with granular and smooth cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

29 March 2022