

CERTIFICATE OF ANALYSIS

Anaplasma phagocytophilum

FA Substrate Slide

Catalog No.:	SLD-IFA-AP
Lot:	P150511-002
Expiration:	16 June 2019

Quality Control Method: Indirect FA using *Anaplasma phagocytophilum* positive control (catalog no. PC-IFA-AP), negative control (catalog no. NC-IFA-AP), anti-equine IgG FITC conjugate (catalog no. CJ-F-EQUG-1ML or 10ML).

Specific Reaction: 2-4+ fluorescence on the positive control and negative on the negative control with no background. There are 0 to 60 infected cells per high-power field.

Other Comments: NA

QC Manager: *Torrey Beeson* Torrey Beeson

QC Release Date: 16 June 2016

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality\COAs and TDS\COA (VMRD)\8-, 10-, 12-well Slides\A. phagocytophila (was E. equi)\SLD-IFA-AP_P150511-002_190616.docx

17 June 2015