

CERTIFICATE OF ANALYSIS

*Anaplasma phagocytophilum*

FA Substrate Slide

Catalog No.:	SLD-IFA-AP
Lot:	P171218-003
Expiration:	16 March 2022

Quality Control Method: Indirect FA using *Anaplasma phagocytophilum* positive control (catalog no. PC-IFA-AP), negative control (catalog no. NC-IFA-AP), anti-equine IgG FITC conjugate (catalog no. CJ-F-EQUG-1ML or 10ML).

Specific Reaction: 2 – 4+ signal with the positive control, no background, and negative with the negative control with trace background. There are 0 - 60 infected cells per high-power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 19 March 2018

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

02 April 2018